

# Rhino Titanium Script Form

Defensive - Low Point (suitable for 0.5-1.5 player)

90900T



**To Order** Fill in the frame dimensions and features required. For variations describe in detail on a separate page. If you have any further questions contact us at the address at the bottom of the page.

**SEND FORM**

If online, click here to send via e-mail

## CLIENT DETAILS

Date \_\_\_\_\_

Name \_\_\_\_\_

Address \_\_\_\_\_

Phone \_\_\_\_\_

Fax / Cell \_\_\_\_\_

Email \_\_\_\_\_

Disability \_\_\_\_\_

Grading \_\_\_\_\_

Weight \_\_\_\_\_

Height \_\_\_\_\_

Shoe Size \_\_\_\_\_

T-Shirt Size \_\_\_\_\_

## DELIVERY DETAILS

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

Required by \_\_\_\_\_

(Delivery is normally 4-6 weeks from order confirmation)

## ADDITIONAL NOTES

(you can type here)

Name \_\_\_\_\_

**FOOTPLATE**

- No Footplate       Alloy Footplate Cover - adjustable **AD11690**       Custom Alloy Footplate Cover **AD11691**  
 Alloy Footplate Cover - fixed **AD11690**       Alloy Footplate Cover - 4 post **AD11692**  
 (Heavy Duty)

**ANTI-TIPS**

- Twin Swivel (wide set) **AD15151**     Space required under Castor \_\_\_\_\_ (max 20mm)  
 Poly Plating on Anti Tip **AD15166**     Anti Tip Protect - Straight Bar \_\_\_\_\_ (height)

**BEARINGS**

- Standard Bearings       Melrose Hybrid Ceramic Bearings **AD20400**

**CASTORS**

- 76 Rollerblade (std)       Other \_\_\_\_\_

**UPHOLSTERY**

- Seat Sling (std)       No Cushion       Hip Strap - webbing **AD16510**  
 Rhino Back Upholstery Adjustable Height only (std) \_\_\_\_\_ (width)     Cushion 1" **AD16201**       Hip Strap - ratchet **AD16550**  
 Rhino Back Upholstery Adjustable Width and Height **AD16050**     Cushion 2" **AD16202**       Foot Strap - webbing **AD16333**  
 Back Bar Cover **AD16420**       Foot Strap - ratchet **AD16350**  
 Leg Strap **AD16330**  
 Chest Strap **AD45210**  
 (measure around your chest and seat back)

- Gloves (style & size) \_\_\_\_\_     Wheel Bag (Black with choice of trim colour) **AD16780** (Choose your trim colour) **Bk R Bl Y Gr Pk**  
 Wheel Bag USA **AD16782**     Classic     Retro     Dirty  
 Custom Wheel Bag **AD16781** (describe below)

**SKIRT GUARDS**

- Fabric Skirt Guards **AD16100**     Alloy Skirt Guards - screw on **AD16130** \_\_\_\_\_ (height)

**ACCESSORIES**

- Ball Holder Mounts **AD20350**  
 Titanium Studs **AD20300**

Seat Width - Rear \_\_\_\_\_

Seat Width - Front \_\_\_\_\_

**A** Seat Length \_\_\_\_\_

**B** Back Width \_\_\_\_\_

**C** Back Splay - Top \_\_\_\_\_

**D** Back Splay - Length \_\_\_\_\_

**E** Back Height \_\_\_\_\_

**F** Back Angle \_\_\_\_\_

**G** Seat to Footplate \_\_\_\_\_

**H** Overall Frame Length  
(centre of camber bar to front of pick bar) \_\_\_\_\_

**I** Wheel Base  
(centre of camber bar to centre of pivot) \_\_\_\_\_

**J** Seat to Floor - Front \_\_\_\_\_

**K** Seat to Floor - Rear \_\_\_\_\_

**L** Balance Point \_\_\_\_\_

**M** Pivot Centres  11.5  12  12.5  13  13.5  14

**N** Footplate Position \_\_\_\_\_

**O** Between Wheels \_\_\_\_\_

**P** Camber \_\_\_\_\_

**WINGS**

- Single Wings **AD20200**
- Double Wings **AD20210T**
- Inverted Wings **AD20220T**

Wing Height \_\_\_\_\_

**PICK BAR**

- Titanium (welded) **AD20132**
- Titanium (bolt on) **AD20131**
- Steel (bolt on) **AD20120**

Width \_\_\_\_\_ Height \_\_\_\_\_

Angle \_\_\_\_\_ Code PB - \_\_\_\_\_

- Q**  Transfer Rail **AD20312**

- R**  Holes for Bindings




- S**  Flexible Backpost Cap **AD20310**  
(only available with fixed back)

- T** Axle Receivers  1/2" (std)

- 5/8" **AD12550**

- Castor Position  **U**  **V**

Name \_\_\_\_\_



Revision 2 Rhino T Def October 2017

Name \_\_\_\_\_

**NO WHEELS** Set the wheelchair up for: \_\_\_\_\_

**WHEELS**

**TYRES**

**24"**

- Spoked Wheels 1/2" axle (std)
- Spinerge SL Wheels 1/2" axle (24 spoke) **AD13410**
- Spinerge SLX Wheels 1/2" axle (24 spoke) **AD13420**
- Spinerge SLX Wheels 5/8" axle (24 spoke) **AD13430**

**24"**

- 24" x 1" High Pressure Smooth (std)
- Kenda Ironcap Blue (std)
- Schwalbe Right Run **AD13740**
- Tufo Tubeless **AD13760**

**25"**

- Spoked Wheels 1/2" axle **AD13060**
- Spinerge SL Wheels 1/2" axle (24 spoke) **AD13410**
- Spinerge SLX Wheels 1/2" axle (24 spoke) **AD13420**
- Spinerge SLX Wheels 5/8" axle (24 spoke) **AD13430**

**25"**

- 25" x 1" High Pressure Smooth (std)
- Kenda Ironcap Blue (std)
- Schwalbe Right Run **AD13740**
- Tufo Tubeless **AD13760**

**26"**

- Spoked Wheels 1/2" axle **AD13060**
- Spinerge SL Wheels 1/2" axle (24 spoke) **AD13410**
- Spinerge SLX Wheels 1/2" axle (24 spoke) **AD13420**
- Spinerge SLX Wheels 5/8" axle (24 spoke) **AD13430**

**26"**

- 26" x 1" High Pressure Smooth (std)
- Kenda Ironcap Blue (std)
- Schwalbe Right Run **AD13740**
- Tufo Tubeless **AD13760**

**PUSHRIMS**

- Anodised Alloy (std)
- Rubber coated (red) **AD13845**

**PUSHRIMS - POSITION**

- Standard
- In close
- Onto rim

**SPOKE GUARDS**

- Nil
- Clear Polycarb Dish Spoke Guard **AD13920**
- Clear Polycarb Attacker Spoke Guard **AD13925**
- Oversized Slotted Black Poly Dish Spoke Guard **AD13935**
- 35mm hole
- Other hole (size) \_\_\_\_\_

**SPOKE GUARDS with PUSHRIMS (ONE PIECE)**

- Nil
- Alloy SG001 (16mm p/r) **AD13950**
- Alloy SG001 (19mm p/r) **AD13950**
- Alloy Poly SG006 (16mm p/r) **AD13955**
- Alloy Poly SG006 (19mm p/r) **AD13955**
- Alloy Poly Oversize SG006os **AD13956**
- Red Rubber Poly SG004 (16mm p/r) **AD13960**
- Red Rubber Poly SG004 (19mm p/r) **AD13960**

**SPARE WHEELS as above**

- Spoked Wheel 24" **AD13507** (quantity) \_\_\_\_\_
- Spoked Wheel 25"/26" **AD13561** (quantity) \_\_\_\_\_
- Spinerge SL Wheel **AD13641** (quantity) \_\_\_\_\_
- Spinerge SLX Wheel 1/2" axle **AD13651** (quantity) \_\_\_\_\_
- Spinerge SLX Wheel 5/8" axle **AD13661** (quantity) \_\_\_\_\_

**TYRES, TUBES and AXLES on SPARE WHEELS**

- Tube fitted **AD13695** (quantity) \_\_\_\_\_
- Kenda Tyre fitted **AD13701** (quantity) \_\_\_\_\_
- Schwalbe Tyre fitted **AD13741** (quantity) \_\_\_\_\_
- Axle 1/2" **AD12508** (quantity) \_\_\_\_\_
- Axle 5/8" **AD12510** (quantity) \_\_\_\_\_

**SPOKE GUARDS for SPARE WHEELS as above**  FIT TO SPARE WHEEL

- Clear Polycarb Dish Spoke Guard **AD13921** (quantity) \_\_\_\_\_
- Clear Polycarb Attacker Spoke Guard **AD13926** (quantity) \_\_\_\_\_
- Oversized Slotted Black Poly Dish Spoke Guard **AD13936** (quantity) \_\_\_\_\_
- Red Rubber Pushrim **AD13846** (quantity) \_\_\_\_\_

DO NOT FIT TO SPARE WHEEL

- Alloy SG001 **AD13951** (quantity) \_\_\_\_\_
- Alloy Poly SG006 **AD13954** (quantity) \_\_\_\_\_
- Alloy Poly Oversize SG006os **AD13957** (quantity) \_\_\_\_\_
- Red Rubber Poly SG004 **AD13961** (quantity) \_\_\_\_\_

**EXTRA AXLES**

- Axles 1/2" **AD12514** \_\_\_\_\_
- Axles 5/8" **AD12516** \_\_\_\_\_

**EXTRA BEARINGS**

- Steel 608 5/16" **AD20402** \_\_\_\_\_
- Steel R8 1/2" **AD20405** \_\_\_\_\_
- Ceramic 608 5/16" **AD20403** \_\_\_\_\_
- Ceramic R8 1/2" **AD20407** \_\_\_\_\_